



ACCREDITED CALIBRATION LABORATORIES

# Pipette Calibration & Preventative Maintenance

[Download Order Form](#)



CATLab offers Calibration, Preventive Maintenance and Repair for all major brands of Single and Multi-Channel Pipettes pursuant to ISO 8655, Certified and Traceable to the National Institute of Standards and Technology.

CATLab technicians use A&D Pipette Testers with evaporation trap to minimize error and improve repeatability. The balance interfaces with Indysoft's validated Pipette Extension Module via RS232 to automatically record test data, convert mass-to-volume and generate your calibration certificate while assuring compliance with 21 CFR Part 11, GMP, GLP and GCP guidelines.

- Accupet
- BioHit
- Brantech
- Corning
- Eppendorf
- Finnipette
- Gilson
- Hamilton
- Labnet
- Matrix
- Mettler Toledo
- Rainin
- Thermo Fisher
- Vistalab
- VWR
- & more

## We Keep it Simple

### Calibration

1. Calibration 10/10% & 10/100%, Record as-found data
2. Calibration, Adjustment, Record as-left data
3. NIST Traceable Certificate w/ AF/AL data + label

### Repair & Calibration

1. Inspect, replace components, grease & reassemble
2. Calibration, Adjustment, record as-left data
3. NIST Traceable Certificate w/ AF/AL data + label

### Preventative Maintenance & Calibration

1. Calibration 10/10% & 10/100%, Record as-found data
2. Remove shaft, clean & inspect components
3. Replace O-Rings & Minor Repair as required
4. Grease & reassemble
5. Calibration, Adjustment, Record as-left data
6. NIST Traceable Certificate w/ AF/AL data + label
  - ✓ Custom calibration points upon request
  - ✓ 6-month warranty on PM & Repair

## CATLab – Accredited Calibration Laboratories

Supporting Process Control Professionals Since 1989

[www.CATLab.net](http://www.CATLab.net)

Contact: [sales@CATLab.net](mailto:sales@CATLab.net) for more info or quote



Request pickup in North Carolina  
[NC.Lab@CATLab.net](mailto:NC.Lab@CATLab.net)

877-4CATLab

Request Pickup in Virginia  
[VA.Lab@CATLab.net](mailto:VA.Lab@CATLab.net)